

IOB Forced-Air Warming System

The Best Choice for Body Temperature Management

Safe, Reliable, Easy to Use



Dimensions: H28cm x D22cm x W22cm

Weight: 5.4 kg

Relative Noise Level: <55 dBA

Operating Temperature:

High: $43^{\circ} \pm 2.0^{\circ}\text{C}$ ($109.4^{\circ} \pm 3.6^{\circ}\text{F}$)

Med: $38^{\circ} \pm 2.0^{\circ}\text{C}$ ($100.4^{\circ} \pm 3.6^{\circ}\text{F}$)

Low: $32^{\circ} \pm 2.0^{\circ}\text{C}$ ($89.6^{\circ} \pm 3.6^{\circ}\text{F}$)


Device Rating:

110-120 VAC, 60 Hz, 12 A

220-240 VAC, 50/60 Hz, 8A

Filter: 0.2 μm High-efficiency Filter, 500-1000 hours

Standards: IEC 60601-1, IEC 60601-1-2, IEC 60601-1-8, IEC 80601-2-35



Medisafe Distribution Inc.

7 - 11771 Horseshoe Way, Richmond, BC V7A 4V4 Canada


Ph: (604) 232-2442 Toll Free: 1-855-246-2442 Fax: (604) 232-2445

Email: sales@medisafecanada.com www.medisafecanada.com


IOB Warming Blanket

Intraoperative Blankets


Upper Body
Model - 003



Lower Body
Model - 002



Torso
Model - 001



Full Body Surgical
Model - 008




Underbody Blankets


Pediatric Underbody
Model - 005




Large Pediatric Underbody
Model - 009




Spinal Underbody
Model - 010



Lithotomy Underbody
Model - 011




Full Access Underbody
Model - 006




Pediatric Blankets

Pediatric Lower Body
Model - 012



Pediatric Long
Model - 014




Pediatric Full Body
Model - 007



Pediatric Underbody
Model - 005




Large Pediatric Underbody
Model - 009




Postoperative Blankets

Adult Full Body
Model - 004



Pediatric Full Body
Model - 007




Other Blankets

Cath Lab
Model - 015



Surgical Access
Model - 016



OTHER STYLES ARE AVAILABLE



Medisafe Distribution Inc.
7 - 11771 Horseshoe Way
Richmond, BC V7A 4V4 Canada

Ph: (604) 232-2442 Toll Free: 1-855-246-2442 Fax: (604) 232-2445
Email: sales@medisafecanada.com www.medisafecanada.com